


## Hydration of two calcium sulfates (plaster and anhydrite III)

### Introduction:

Four forms of calcium sulfate are known : dihydrate  $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$  (gypsum), hemihydrate  $\text{CaSO}_4 \cdot \frac{1}{2}\text{H}_2\text{O}$  (plaster) and two anhydrites  $\text{CaSO}_4$  : natural and insoluble anhydrite and soluble anhydrite III (obtained by heating plaster to  $200^\circ\text{C}$ ).

Plaster and anhydrite III, which are soluble in water, are used in the building industry, because of their setting properties.

Calorimetry is an interesting method for studying the setting of calcium sulfates.



### Experimental:

- Samples :
- 1) Plaster (625 mg)
  - 2) Anhydrite III (625 mg)
  - 3) Water (500 mg)

Crucible : Reversing mixing cell

Heating mode : Isothermal  $28^\circ\text{C}$

### Conclusion:

At initial time, sample and water are separated by a lid in the mixing cell.

The mixing of the two compounds is carried out by reversing the calorimeter.

Two stages are seen on the recorded thermograms :

- . dissolution of the sample in water
- . setting of the mixture (hydration)

The setting can be characterised by the half crystallization time (at the top of the hydration peak).

In this case the setting of plaster ( $t_{1/2} = 25 \text{ min}$ ) is faster than that of anhydrite ( $t_{1/2} = 50 \text{ min}$ ).

### Instrument:

C80

Ambient up to  $300^\circ\text{C}$



[www.setaram.com](http://www.setaram.com) – sales @ [setaram.com](mailto:sales@setaram.com)